

ergo-Logistcis

ergo-Logistcis.com

Customer Service:

Joan Hinchman

1-800-341-7048

Fax: 1-800-704-2815

CustomerService@ergo-Logistics.com



Technical Support: 509-747-5027

email: Support@ergo-Logistcis.com

Website:

<http://www.ergo-logistics.com/>

Countertop Washer Decontam **PASS-THRU** or Single Door

STEELCO Washers DS-50, DS-50-DRS, and DS-50-DRS-2

STEELCO Washers DS-50, DS-50-DRS, and DS-50-DRS-2

External dimensions 632x562x535mm [24.9"x22.12"x21.06"] (with door closed)

External dimensions 632x852x535mm [24.9"x33.54"x21.06"] (with open door)

Chamber usable Size Wide 19.69" x 15.75" x Deep 9.65 High

STEELCO WASHER DS-50

Countertop Tabletop Washer

STEELCO WASHER DS-50-DRS

Countertop Tabletop Washer with **DRYER**

STEELCO WASHER DS-50-2 DRS **2 door Pass Through**

the **only** Countertop Tabletop Pass Through Washer with **DRYER**

208 volt or 220 volt, 60Hz, UL Listed

STEELCO DS50, DS50DRS, and DS50DRS-2 RACKS

Measurements - W x D x H

Incert holds 35 Instruments vertically 395 x 125 x 120mm [15.55"x4.92"x4.72"] **C03**

Incert holds 8 Cassettes 36mm [1.4"] 390 x 180 x 110mm [15.35"x7.08"x4.33"] **C04**

Primary Rack 1 Level with Sides - 470 x 390 x 85mm [18.5"x15.35"x3.34"] **C06**

Cleaning Basket w Lid S/S perforated 100 x 100 x 40mm [3.93x3.93"x1.57"] **C28**

Insert for Instrument Containers - for holding 3 Cassettes 11 x 8 x 1.5 vertically **C90**

Insert for Instrument Containers - for holding 4 Cassettes 11 x 8 x 1.5 vertically **C91**

Tray for instruments **C92**



C03



C04



C06



Glass rack 18 position	C83
Rack with 4 pipes & hand piece connections	C84
Rack for dental instruments with 8 places	C85
Insert for instrument containers with 3 positions, Cassettes 280x200x40h mm	C90
Insert for instrument containers with 4 positions, Cassettes 205x120x40h mm	C91
Big tray for instruments, Delivery cassettes, 280x200x40h mm	C92
Rack for hand pieces with 8 places DENTAL	C640
Dental Handpiece support - lower part - brochure page 15	C080009
Dental Handpiece support - upper part - brochure page 15	C080010
Dental Handpiece support - stainless steel filter - brochure page 15	C080006
Dental Handpiece support - adapter Ø 20mm - brochure page 15	C080005
Dental Handpiece support - adapter Ø 15mm - brochure page 15	C080001

C28

STEELCO DS50, DS50DRS, and DS50DRS-2

Pump to Drain. Does not require a dedicated floor drain. Can pump into Sink Drain Plumbing or Wall Drain.

STEELCO DS50, DS50DRS, and DS50DRS-2 Dimensions: (W x D x H)

External dimensions 632x562x535mm
[24.9"x22.12"x21.06"] (with door closed)

External dimensions 632x852x535mm
[24.9"x33.54"x21.06"] (with open door)

Max rack dimensions 400x500x245mm
[15.75"x19.69"x9.65"]

Max. basket capacity 49 lt. [13gal.]

Electrical connection 1x208V/60Hz

Electrical tank heating 3 kW

Washing circuit with one pump for a total of 110W

Total power 3.5 kW

Water connections: one

STEELCO DS50, DS50DRS, and DS50DRS-2 Description:

washing chamber entirely made of stainless steel AISI 316L

smooth chamber design to eliminate dead legs for potential bacteria grow

self cleaning chamber with rounded corners

water filtering system for excellent washing results

OPTIONAL Pass through version with two manual bottom hinged doors

Frontal basket loading

Automatic Rack coupling system for injector cleaning and drying to guarantee highest quality

cleaning and efficient drying

cleaning and drying results

- Full drainage of re-circulating pump
- Thermal decontamination possibilities up to 93°C [200°F]
- PT1000 temperature probe
- Equipped with 2 dosing pumps for chemicals, microprocessor controlled.
- Level control for chemical jugs with minimum level alarm
- LED display as user interface
- 20 Pre-set programs

Total control of cycle parameters:

time, temperature, chemical dosing, water quantity,

drying

- Exhaust air condenser
- Equipped with hot air dryer

not included in price:

- > start-up and training
- > installation
- > transport & insurance

STEELCO DS 500 DRS Warranty: 1 year parts

not included in price:

- > start-up and training
- > installation
- > transport & insurance