

all-in-ONE™ Enzymatic Enzyme Detergent Surgical Instrument Cleaner with Conditioners

Highly concentrated, fast, and cost effective presoak and cleaner.
Mallard Medical website: www.MallardMedicalSurgicalInstrumentCare.com

Image: replacing multiple products, lowering costs, and improving performance with the all-in-ONE™ Surgical Instrument Cleaner AND Conditioners. The all-in-ONE™ rapidly breaks down and removes all forms of proteinaceous bioburden, cleans and conditions the surface, enhances the “passive layer” protection of surgical stainless steel, lubricates moving parts, and then rinses clean. The all-in-ONE™ highly concentrated Surgical Instrument Cleaners AND Conditioners will deliver **512 gallons** (1938 liters) of Surgical Instrument Cleaning AND Conditioning power. The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and **STAINS**. The **Neutral pH** all-in-ONE™ is non-irritating to eyes, skin, and is **100% biodegradable**.

The all-in-ONE™ easy FOAM-it™ lays down a thick blanket of long lasting foam to hydrate and **remove** all forms of proteinaceous debris, hard water deposits, and encrusted **STAINS**, for fast acting cleaning action. When water is added the ‘foam’ dissipates providing a **clear ‘free-rinsing’** cleaning solution

all-in-ONE™ Surgical Instrument Cleaner with Conditioners

ONE product for Cleaning and Conditioning surgical instruments and scopes.

ONE application for thoroughly Cleaning and Conditioning.

Apply it **ONCE...**

It's all you really need for Cleaning and Conditioning Surgical Instruments and Scopes.

Replace Cleaning and Conditioning products for: removing all forms of proteinaceous bioburden, pre-cleaning and soaking to remove **STAINS**, manual cleaning, Ultrasonic cleaning, and automated Surgical Instrument Washing.



easy FOAM-it™
Surgical Instrument
Cleaner and Conditioners



all-in-ONE™
Surgical Instrument
Cleaner with Conditioners

The all-in-ONE™ Surgical Instrument Cleaners deliver a unique formulation of detergents and multi-tiered "high-level" enzymes (Protease, Amylase, Lipase, Carbohydrase) within a **Neutral pH** surfactant chemical complex. The all-in-ONE™ cleaners will: break down bioburden and clean the surface, remove **STAINS**, enhance the "passive layer" protection of surgical stainless steel, and lubricate moving parts. The **highly concentrated all-in-ONE™** cleaners are aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment. Compare the cleaning performance of the all-in-ONE™ to the lesser effective 'enzyme' detergents, i.e. the single and dual enzyme (one and two enzyme) detergents STERIS Prolystica Enzymatic Presoak and Cleaner and the MetriZyme Dual Enzymatic Detergent. You will discover the superior cleaning performance of the all-in-ONE™ Surgical Instruments Cleaner with Conditioners.

UPS-FS-32 all-in-ONE™ easy FOAM-it™ 12/case
Ergo UPS1-1 all-in-ONE™: 4 gallons/case
Ergo UPS1-5 all-in-ONE™: 1 ea. 5 gallon
Ergo UPS1-15 all-in-ONE™: 1 ea. 15 gallon
Ergo UPS1-30 all-in-ONE™: 1 ea. 30 gallon
Ergo UPS1-55 all-in-ONE™: 1 ea. 55 gallon

Dilution Rate: .25 ounce to .5 ounce per gallon depending on the level of bio-burden and **STAINS**.

Mallard Medical: 1-800-441-1907 FAX: (215) 297-5920 email: MallardMedical@comcast.net
website: www.MallardMedicalSurgicalInstrumentCare.com
Technical Support: Support@ergo-Logistics.com 1-509-747-5027 Website: ergo-Logistics.com