



Customer Service: 1-509-747-5027, email: [CustomerService@ergo-Logistics.com](mailto:CustomerService@ergo-Logistics.com)  
Technical Support: John Temple, 1-509-981-1553, email: [Support@ergo-Logistics.com](mailto:Support@ergo-Logistics.com)

## Bedside Care Kit, Endoscope Cleaning Sponges, Sponges Surgical instrument Cleaning

### Surgical Instrument Cleaning Sponges

#### Endoscope Cleaning Sponges

#### Sponge Cleaning Kit for Medical Devices and Equipments

The D-Zyme Sponge retains more enzyme detergent than other sponges and paper based wipes.

The all-in-ONE 4 Enzyme Detergent Cleaner is over 4 times more effective than other enzyme solutions, delivering the 4 Enzymes

(Protease, Amylase, Carbohydrase, and Lipase Enzymes) needed to breakdown all forms of proteinaceous bioburden

PLUS the surface cleaning Detergents needed for rapid and effective cleaning.

Non-irritating, 100 % Biodegradable, Neutral pH, and Phosphate free.

### Surgical Instrument Cleaning Sponges

D-Zyme Sponge, Order Number: ESW105200

Flat, Multi-Use, Sponge Loaded with highly concentrated all-in-ONE Enzyme Detergents.

Individual Packaged - 200 ea. Per Case (4 boxes of 50)

Contents: Sponge Pre-soaked with Enzyme Detergent

Sponge Size: 4" Wide x 5" Long x 1/2" High

Order Number: ESW105200



### Endoscope Cleaning Sponges

D-Zyme Sponge, Order Number: ESW104100

Tubular, Multi-Use, Sponge, Loaded with highly concentrated all-in-ONE Enzyme Detergents.

Individual Packaged - 100 ea. Per Case (4 boxes of 25)

Contents: Sponge Pre-soaked with Enzyme Detergent

Sponge Size: 4.5" Long x 2.4" Diameter

Order Number: ESW104100

### Sponge Cleaning Kit for Medical Devices and Equipments

Appli-Kit, Order Number: APT215035

Pouch contains Sponge and 215 ml of the highly concentrated all-in-ONE 4 Enzyme Detergents.

Pouch has a wide base pouch for stability and easy fluid transfer.

Individual Packaged - 35 ea. Per Case

Contents: Pouch with Sponge and 215 ml of Enzyme Detergent

Sponge Size: 4.5" Long x 2.4" Diameter

Order Number: APT215035

